

Release notes
for Permobil R-net PC-programmer ver 6.3.0
Dealer access level

2019-09-26/PH

Updates, changes and new functions in this release:

(R-net Application v33, Parameter info ver 1.1248, Executable ver 6.3.8)

- Nodes now display the correct name under System Logs.

Release notes
for Permobil R-net PC-programmer ver 6.1.5
Dealer access level

2018-12-21/JWL

Updates, changes and new functions in this release:

(R-net Application v33, Parameter info ver 1.1240, Executable ver 6.1.5)

- Support for Omni2
- Support for AIM

Release notes
for Permobil R-net PC-programmer ver 6.1.3
Dealer access level

2016-10-20/JÅF

Updates, changes and new functions in this release:

(R-net Application v33, Parameter info ver 1.1165, Executable ver 6.1.3)

- CJSM2 BT and IR support
- Softkeys introduced including icons and can be set per mode.
- Profile/Mode jack can be set to stereo or mono (for assignable buttons), can be used for Bluetooth or as Assignable Buttons.

Release notes
for Permobil R-net PC-programmer ver 6.0.0
Dealer access level

2016-07-06/JÅF

Updates, changes and new functions in this release:

(R-net Application v32, Parameter info v 1.113912, Executable v 6.1.0)

- Assignable Buttons functionality introduced
- Support for CJSM2/PJSM

Release notes
for Permobil R-net PC-programmer ver 5.5.1
Dealer access level

2014-12-09/JÅF

Updates, changes and new functions in this release:

(R-net Application v32, Parameter info v 1.1031, Executable v 5.3.1)

- Support for iDevice BlueTooth Module
- ICS Limits updated for Corpus AP and Corpus VS

Release notes
for Permobil R-net PC-programmer ver 5.2.1
Dealer access level

2014-04-31/JÅF

Updates, changes and new functions in this release:

(Application v 5.2.1, Parameter info v 1.958, Executable v 5.2.5)

- Support for 2nd BlueTooth Mouse Module
- IR Configurator Library updated to V6

Release notes
for Permobil R-net PC-programmer ver 5.1.0
Dealer access level

2013-04-26/JÅF

Updates, changes and new functions in this release:

(Application v 5.1.0, Parameter info v 1.916, Executable v 5.1.4.12)

- Permobeat parameters removed
- The parameter *Displays* added
- Input device Octopus is now called *Magic Drive Touch*
- When loading a file with higher level (for example OEM-level), there are now the options to View the file or Download the file. If the view option is selected, the file is opened in a read-only mode where editing and Downloading of the file is prohibited. If the option to download the file is selected, the file is downloaded at its current access level.
- IR Configurator Library raised from V3 to V4

Release notes
for Permobil R-net PC-programmer ver 5.0.0
Dealer access level

2012-08-23/JÅF

Updates, changes and new functions in this release:
(Application v 5.0.0, Parameter info v 1.899, Executable v 5.1.3)

- New BT Mouse Parameters
- New IR Config. Library v3
- Creates a shortcut to Start – Permobil – R-net
- Some default ICS parameters settings corrected.

Release notes
for Permobil R-net PC-programmer ver 4.2.0
Dealer access level

2010-11-29/JÅF

Updates, changes and new functions in this release:
(Application v 4.2.0, Parameter info v 1.725, Executable v 4.1.7)

- ICS Parameters added under Seating - ICS
- New Personality files available under ‘Controller’ – ‘Set Controller Personality’
- Clear System Log now possible
- Error Codes displayed correctly (for ICS, ESP, Copilot, HMC-joysticks and devices)

Release notes
for Permobil R-net PC-programmer ver 4.0.2
Dealer access level

2010-04-22/JÅF

Updates, changes and new functions in this release:
(Application v 4.02, Parameter info v 1.60121, PGDT version name 4.04)

- IR-configurator included (for OMNI)
- OMNI edit bug fixed
- Works better on Windows 7
- Some new parameters can now be seen and edited

Release notes

for Permobil R-net PC-programmer ver 3.2.9
Dealer access level

2009-09-14/JÅF

Updates, changes and new functions in this release:

- 'Saved as OEM' files will update Preset data in controller when executing 'write to controller' command
- Improved communication PC – dongle
- Help information structure corrected

Release notes for Permobil R-net PC-programmer ver 3.2.3 Dealer access level

2009-04-15/JÅF

Updates, changes and new functions in this release:

- A parameter file saved with an OEM-programmer (file:OEM) is now directly programmed to controller. The file is no longer opened and settings can't be seen on the PC-screen. An informative text appears on the screen when trying to open an OEM-file with a Dealer programmer.
- New input devices that can be connected to the Rnet bus are now added and accessible, for example Mouse controller, MiniJoystick, iPod and Remote Button Interface.
- Complete back-up of all parameters can now be done with the Dealer programmer via the 'File' – 'Save As OEM'. ('File' – 'Save' only saves Dealer accessible parameters)
- It is now possible to update the OBP-structure. Selectable is Standard ('PGDT Generic') or 'Permobil ESP'. With 'Permobil ESP' becomes the ESP speed parameters accessible via OPB (substituting the normal speed parameters).
Click on 'Controller' – 'Set Controller Personality', select OBP set-up and click OK.
- 'Allow grab' parameter is now possible to change for all profiles with Dealer version.
- Both .r-net and .rnss files are shown in 'File' – 'Open' windows as default.
- Output signal from ESP is possible to view in a similar way as the joystick output via 'Tools' – 'System Tests' – 'Intermediate XY'.
- When installing this version of Permobil PC-Programmer will not the old version be uninstalled. (Should be done manually to avoid using the old one). That means that the Permobil Dealer PC-programmer can exist together with any other Dealer PC-programmer (for example PGDT's generic version) on the same PC without interference. Useful for dealers and others working with different brands of wheelchairs.
- Some more bugs are fixed, such as bugs related to input type and other profile management settings and handling (e.g. profiles effected by subsetfiles when should not).
- Complete installation file for this version is now stored on the Permobil webb-site. (Not only exe-file and rnd-file etceteras as before). Installation file is >20Mb.

Svenska

**Release notes
för Permobil R-net PC-programmer ver 4.2.0
Dealer access level**

2010-11-29/JÅF

(Sv) Uppdateringar, ändringar och nya funktioner i denna release:
(Application v 4.2.0, Parameter info v 1.725, Executable v 4.1.7)

- ICS Parameters tillagt under Seating - ICS
- Nya Personality- filer under 'Controller' – 'Set Controller Personality'
- Clear System Log är nu möjligt.
- Felkoder visas nu korrekt (för ICS, ESP, Copilot, HMC-joysticks och enheter)

**Release notes
för Permobil R-net PC-programmer ver 4.0.2
Dealer access level**

2010-04-22/JÅF

(Sv) Uppdateringar, ändringar och nya funktioner i denna release:
(Application v 4.02, Parameter info v 1.60121, PGDT version name 4.04)

- IR-configurator inkluderad (för OMNI)
- OMNI edit bug fixad
- Fungerar bättre med Windows 7
- Några nya parametrar kan visas och ändras

**Release notes
för Permobil R-net PC-programmerare ver 3.2.9
Dealer-version**

2009-09-14/JÅF

(Sv) Uppdateringar, ändringar och nya funktioner i denna release:

- Filer som sparas med 'Saved as OEM' kommer nu att uppdatera Preset Data när de skrivs till till systemet med 'Write to controller'.
- Förbättrad kommunikation PC – dongle
- Help information rättad.

**Release notes
för Permobil R-net PC-programmerare ver 3.2.3**

Dealer-version

2009-04-15/JÅF

(Sv) Uppdateringar, ändringar och nya funktioner i denna release:

- En parameterfil som sparats med en OEM-programmerare (file:OEM) skrivs nu direkt till kontrollern. Filen öppnas inte och settings visas inte på PC-skärmen. En informationstext om detta visas när man försöker öppna en OEM-fil med Dealer-programmerare.
- Nya enheter som kan anslutas till Rnet bussen har tillkommit och kan programmeras. Exempel på enheter är Mouse controller, MiniJoystick, iPod och Remote Button Interface.
- En fullständig back-up av alla parametrar kan nu göras med Dealer-programmeraren via 'File' – 'Save As OEM'. ('File' – 'Save' sparar endast Dealer-parametrar)
- Det är nu möjligt att uppdatera OBP-strukturen. Val mellan Standard ('PGDT Generic') eller 'Permobil ESP' kan göras. Med 'Permobil ESP' blir 'speedparametrarna för ESP åtkomliga via OBP (ersätter de vanliga speedparametrarna). Klicka på 'Controller' – 'Set Controller Personality', välj OBP set-up och klicka OK.
- 'Allow grab' parametern kan nu ändras i alla profiler med Dealer-versionen.
- Både .r-net and .rnss filer visas nu i 'File' – 'Open' fönstret som standard.
- Utsignalen från ESPn är nu möjlig att studera på samma sätt som joystickens utsignal via 'Tools' – 'System Tests' – 'Intermediate XY'.
- Vid installation av denna version av Permobil PC-Programmerare kommer inte den gamla att bli avinstallerad. (Bör göras manuellt för att undvika att använda den gamla). Det innebär också att Permobil Dealer PC-programmerare kan existera tillsammans med annan Dealer PC-programmerare (ex PGDT's Generic) på samma PC utan att påverka varandra. Detta är speciellt bra för tekniker/terapeuter som hanterar olika fabrikat av stolar.
- Ytterligare buggar är åtgärdade, såsom ex buggar rörande 'input type' och andra 'Profile Management' parametrar (ex profiler som påverkats av subsetfiler när de inte borde)
- Komplet installationstil för denna version finns nu tillgänglig på Permobil's hemsida. (Inte bara .exe-filen och .rnd-filen etc som förut). Installationsfilen är >20Mb